



sanotracin RTU™

Disinfectant • Bactericide • Sporicide • Tuberculocide • Virucide • Fungicide

Kills 99.9999% *Clostridium difficile* spores in 3 minutes
Disinfects against Acinetobacter baumannii in 2 minutes
Disinfects against MRSA in 2 minutes
Disinfects against VRE in 2 minutes

- Ready-to-use
- Cleans, disinfects, and deodorizes
- Disinfects hard, nonporous surfaces in Hospitals and Health Care Facilities
- Prevents mold, mildew stains and odors
- Does NOT contain Bleach or Alcohol**
- Will not harm LCD screens

ACTIVE INGREDIENTS:

Hydrogen Peroxide.....	4.40%
Peroxyacetic Acid.....	0.23%
Other Ingredients.....	95.37%
Total.....	100.00%

PN 81 00000008F

KEEP OUT OF REACH OF CHILDREN

WARNING ADVERTENCIA

Si usted no entiende la etiqueta, busque a alguien para que se la explique a usted en detalle.
If you do not understand the label, find someone to explain it to you in detail.

See back panel for additional precautionary and first aid statements.

NET CONTENTS: 1 U.S. Quart (.946 Liters)

Do Not Remove
Product Labeling



sanotracin RTU™



PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS

WARNING. Causes substantial but temporary eye injury. Do not get in eyes or on clothing. Wear protective eyewear such as goggles, face shield, or safety glasses and gloves. Harmful if swallowed. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove and wash contaminated clothing before reuse. **PHYSICAL AND CHEMICAL HAZARDS:** This product contains an oxidizing agent. Avoid contact with other sanitizers, cleaners or organic substances. **ENVIRONMENTAL HAZARDS:** This pesticide is toxic to birds, fish and aquatic invertebrates. Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans, or other waters unless in accordance with the requirements of a National Pollutant Discharge Elimination System (NPDES) permit and the permitting authority has been notified in writing prior to discharge. Do not discharge effluent containing this product to sewer systems without appropriate dilution. **FIRST AID: IF IN EYES:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice. **IF SWALLOWED:** Call a poison control center or doctor immediately for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to by a poison control center or doctor. Do not give anything to an unconscious person. Have the product container or label with you when calling a poison control center or doctor or going for treatment. **For emergency medical information, call 1-888-779-1339.**

Sanotracin RTU™ solution will not harm aluminum, LDPE, HDPE, vinyl, painted surfaces, polycarbonate, polypropylene, polyurethane varnish, PVC, silicone rubber, stainless steel, medical tubing, vinyl rubber, acrylic, brass, LCD screens, copper, Corian®, laminate flooring and countertops, nickel, glazed porcelain, glass, and glazed tile.

Sanotracin RTU™ is a trademark of BioMed Protect, LLC, St. Louis, MO, USA. Patents pending.

EPA Reg. No. 88089-4-12120
EPA Est. No. 88089-MO-1^A, 87720-MO-1^A
67811-MO-1^A, 4313-TX-1^P

Superscript is first letter of lot number.
See Lot No. on Bottle.

For Emergency Medical Information:

1-888-779-1339

For Product Information: 1-800-323-2251

www.triple-s.com

PRODUCT NO.

18001



PN 810000008_2

Manufactured for:



2 Executive Park Dr.
BILLERICA, MA 01862

STORAGE AND DISPOSAL Do not contaminate water, food, or feed by storage or disposal. **STORAGE:** Store in original containers in a cool, dry area, away from direct sunlight. Avoid prolonged storage temperature above 40°C, 104°F. This may cause increased degradation of the product, which will decrease product effectiveness. In case of spill, flood area with large quantities of water. **DISPOSAL:** Wastes resulting from the use of this product must be disposed of onsite or at an approved waste disposal facility. **CONTAINER DISPOSAL: Non-refillable container. Do not reuse or refill this container.** If empty, wrap container and put in trash or offer for recycling, if available. If partly filled, call your local solid waste agency for disposal instructions.

WARRANTY This material is safe and effective under the conditions of application as specified in the Directions for Use. To the fullest extent permitted by law, the manufacturer shall not be liable for improper application or use of this product beyond those purposes referred to in the Directions for Use. NO OTHER EXPRESS OR IMPLIED WARRANTY OF FITNESS OR MERCHANTABILITY IS MADE.

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling. Sanotracin RTU™ microbicide can be used for disinfection, mold and mildew control. This product must be used only as directed in this label. This product is a "ready-to-use" liquid, comprised of 4.4% concentration of hydrogen peroxide and 0.23% concentration of peroxyacetic acid, and it should be used at full strength and should not be diluted before use. **Before using this product read the entire label.**

Sanotracin RTU™ can be applied to hard, non-porous surfaces such as counter-tops, tables, furniture, cabinets, appliances, kitchen and bathroom fixtures, furnishings, hardware, floors, walls, equipment, tabletops, benchtops, drip pans, troughs, desks, and locker room benches. For heavily soiled areas, pre-cleaning is required. Apply Sanotracin RTU™ microbicide to the contaminated surface with a spray applicator device at an application rate greater than 50 ml/m² (approximately 1.5 fluid oz. of the "ready-to-use" Sanotracin RTU™ liquid per square yard), or by other means to ensure that the surface is wet (for example, by wiping with a clean cloth or sponge that is wet with Sanotracin RTU™ liquid), or by soaking the contaminated surface and allowing it to remain on the surface for the indicated contact time. Spray or soak the surface with a sufficient amount of Sanotracin RTU™ solution to remain on the surface for the indicated contact time before evaporation. Air dry or remove with a dry clean cloth or wipe.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

DISINFECTION Use Sanotracin RTU™ microbicide as a disinfectant in hospitals, surgical and obstetrical suites, housekeeping services, physical therapy departments, nursing services, autopsy facilities, nursing homes, health-care facilities, schools, colleges, day-care centers, child-care facilities, veterinary clinics, animal life science laboratories, industrial facilities, dining areas (non-food contact surfaces), office buildings, recreational facilities, retail and wholesale establishments, and in transport vehicles including ships, planes, automobiles, trucks, trains, ambulances, and EMS.

***Combination Disinfection and Cleaning:** Sanotracin RTU™ is effective against *Staphylococcus aureus*, *Salmonella enterica*, *Enterobacter aerogenes*, *Salmonella typhimurium*, *Listeria monocytogenes*, *Vibrio cholerae*, *Pseudomonas aeruginosa*, and other bacteria listed in the Chart, on hard, non-porous surfaces. For heavily soiled areas a pre-cleaning step is required. Apply Sanotracin RTU™ with a spray applicator device or other means that ensures that the surface is thoroughly wet. Allow surfaces to remain wet for **2 minutes**; then, wipe with a dry, clean cloth. If the Sanotracin RTU™ becomes soiled (i.e., is not visibly clear), or diluted, replace and use fresh Sanotracin RTU™. **Virucidal:** Sanotracin RTU™ is effective against *Influenza virus Type A (H3N2)*, *Influenza virus Type A (H1N1)*, *Norovirus*, *Poliovirus*, *Reovirus*, *Rhinovirus*, *Human coronavirus*, and other specific viruses listed in the Chart, on hard, non-porous, inanimate surfaces with a **2 minute** contact time in the presence of a 5% organic load. For heavily soiled areas a pre-cleaning step is required. Apply Sanotracin RTU™ with a spray applicator device or other means that ensures that the surface is thoroughly wet. Allow surfaces to remain wet for **2 minutes**; then, wipe with a dry clean cloth. **Fungicidal:** Sanotracin RTU™ can be used on hard non-porous surfaces such as shower room floors, shower stalls, locker room benches, and hard non-porous bathmats (plastic or rubber). At full strength (4.4% concentration of hydrogen peroxide and 0.23% concentration of peroxyacetic acid), this product is effective against *Trichophyton mentagrophytes* (athlete's foot fungus), *Aspergillus niger* (mold), and *Candida albicans* (pathogenic yeast) in the presence of a 5% organic load with a **1 minute** contact time. For heavily soiled areas a pre-cleaning step is required.

Apply Sanotracin RTU™ with a spray applicator device or other means that ensures that the surface is thoroughly wet. Allow surfaces to remain wet for **1 minute**; then, wipe with a dry clean cloth.

Sporicidal Disinfectant: Sanotracin RTU™ is effective against *Clostridium difficile* spores on hard, non-porous inanimate surfaces when used at full strength (4.4% concentration of hydrogen peroxide and 0.23% concentration of peroxyacetic acid), with a **3 minute contact** time. Apply Sanotracin RTU™ with a spray applicator device or other means that ensures that the surface is thoroughly wet. Allow surfaces to remain wet for **3 minutes**; then, wipe with a dry clean cloth. **Personal Protection:** Wear appropriate barrier protection such as gloves, gowns, masks and eye covering. **Cleaning Procedure:** Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with a clean cloth, mop, and/or sponge saturated with the disinfectant product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of spores. Restrooms are to be cleaned last. Do not reuse soiled cloths.

Infectious Materials Disposal: Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal. **Specific Instructions for Blood-borne Pathogens:** KILLS HIV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in healthcare settings and other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/object likely to be soiled with blood/body fluids can be associated with the potential for transmission of *Human immunodeficiency virus* Type 1 (HIV-1) associated with AIDS.

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATING AGAINST HIV-1 (AIDS Virus) OF SURFACE/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.

Personal Protection: Disposable latex or vinyl gloves, gowns, face masks, or eye coverings as appropriate, must be worn during all cleaning of body fluids, blood, and decontamination. **Cleaning Procedure:** Blood/body fluids must be thoroughly cleaned from surfaces/objects before application of this product. **Contact Time:** Allow surface to remain wet for 2 minutes to kill HIV-1. Use contact times listed for all other organisms claimed. **Disposal of Infectious Material:** Blood/body fluids should be autoclaved and disposed of according to Federal, State, and local regulations for

infectious waste disposal. **Tuberculocidal:** Sanotracin RTU™ microbicide is effective against *Mycobacterium bovis* on hard, non-porous inanimate surfaces with a **5 minute** contact time in the presence of a 5% organic load. For heavily soiled areas a pre-cleaning step is required. Apply Sanotracin RTU™ with a spray applicator device or other means that ensures that the surface is thoroughly wet. Allow surfaces to remain wet for **5 minutes**; then, wipe with a dry clean cloth.

Effective against a wide range of bacteria, viruses and fungi including:

Bactericidal (2 minutes): *Enterobacter aerogenes* (ATCC 13048); *Listeria monocytogenes* (ATCC 19111); *Pseudomonas aeruginosa* (ATCC 15442); *Salmonella typhimurium* (ATCC 13311); *Staphylococcus aureus* (ATCC 6538); *Vibrio cholerae* (ATCC14035); *Salmonella enterica* (ATCC 10708); *Acinetobacter baumannii* (ATCC 19606); *Campylobacter jejuni* (ATCC 29428); *E. coli* O157:H7 (ATCC 35150); *E. coli* ESBL (ATCC BAA-196); *Enterococcus faecalis* - Vancomycin Resistant (VRE) (ATCC 51299); *Enterococcus hirae* (ATCC 10541); *Haemophilus influenzae* (ATCC 19418); *Klebsiella pneumoniae* (ATCC 4352); *Legionella pneumophila* (ATCC 33153); *Proteus vulgaris* (ATCC 27973); *Serratia marcescens* (ATCC 13880); *Shigella sonnei* (ATCC 11060); *Methicillin Resistant Staphylococcus aureus* (MRSA) (ATCC 33592); *Community Acquired - Methicillin Resistant Staphylococcus aureus* (CA-MRSA) (Genotype 300), *Streptococcus pneumoniae* - Penicillin Resistant (PRSP) (ATCC 700671), *Streptococcus pyogenes* (ATCC 19615)

Tuberculocidal (5 minutes): *Mycobacterium bovis*

Virucidal (2 minutes): *Influenza virus* Type A (H3N2) (ATCC VR-2072); *Human coronavirus* (ATCC VR-740); *Adenovirus*, Type 2 (ATCC-846); *Avian influenza virus* H5N1 (NIBRG-14); *Herpes simplex virus*, Type 1 (ATCC VR 260); *Herpes simplex virus*, Type 2 (ATCC VR 734); *Human immunodeficiency virus* Type 1 (HIV-1, Zeptomatrix Corp.), *Influenza virus*, Type A (ATCC VR 544); *Influenza A* (H1N1) (ATCC VR 1469, Strain A/PR/8/34); *Influenza virus*, Type B (ATCC VR 823); *Norovirus* (F. Calicivirus) (ATCC VR 782); *Poliovirus* Type 1 (ATCC VR 1562); *Reovirus* (ATCC VR 232); *Respiratory syncytial virus* (RSV) (ATCC VR-26); *Rhinovirus* (ATCC VR-338); *Rotavirus* (Strain WA)

Sporicidal Activity of Disinfectant (3 minutes): *Clostridium difficile* spores (ATCC 43498)

Fungicidal/Mold Killing (1 minute hard surfaces): *Aspergillus niger* (ATCC 16404); *Candida albicans* (ATCC 10231); *Trichophyton mentagrophytes* (athlete's foot fungus) (ATCC 9533)

Mold Killing (10 minutes) Porous Surfaces: *Aspergillus niger*

PRESS TO RESEAL